

A new way. A better way.

Dental Imaging Software – Reimagined





Capture, store, and access your favorite smiles to and from the cloud in real-time.

You bet your iris! We've simplified things by pre-setting image manipulations, capturing digital images, and securely storing them in the cloud in real-time.

Robust imaging features can be easily added to Curve Hero, eliminating external bridges or cumbersome third-party integrations.

Key Imaging Features



Ease of use

- By overlaying FMX images over patients' restorative charts, you get a clear understanding of what you're examining
- Easily control the opacity with a sliding bar to view the image with or without the chart
- All images are integrated with the patient record, so you can access all the information you need in one place



Ease of Accessibility

- Easily access images from anywhere via a Chrome internet browser
- All images are safely backed up to the cloud in real time
- Data storage and cybersecurity are included in your monthly Curve Dental subscription



Image manipulations

- Save time by recording and saving a series of preferred manipulations and filters
- Apply a saved set of adjustments to any images with a click of the mouse
- Easy adjustments for contrast, brightness, gamma, rotate, sharpen, and invert



Measurement tools

- Measure root depth angles, signifying impingement, areas, and distances for abnormalities and pockets
- Trace paths through the anatomy to establish the right treatment approach
- Save all measurements, notes, and images directly to the patient record













Imaging Compatibility

Curve understands the complexity of your infrastructure and the hardware investments you've made. We work hard to deliver you with flexibility by continually enhancing Curve Hero's compatibility with leading brands.

Digital X-ray Sensors

Carestream 6200

Cyber Medical Imaging XDR

DEXIS Titanium by KaVo

DEXIS Platinum

Gendex GXS-700 (blue sensor)

Image Works (EVA)

Jazz Imaging Solo

KaVo IXS

And more

Masterlink EDLENi Apex

Planmeca ProSensor USB**

Schick 33

Schick CDR with HS Remote

(white box only)

Schick Elite with Elite Remote or HS

Remote

Schick WiFi

Vatech EzSensor

Vatech EzSensor P

Video Dental Concept QuickRay

2D Digital Pans

CS 8100 Digital Panoramic System

CS 9300 System

Planmeca ProMax 2D S2

Planmeca ProMax 2D S3

Planmeca ProOne

Sirona ORTHOSPHOS XG 3

Sirona ORTHOSPHOS XG 5

Sirona ORTHOSPHOS XG 3DReady

Sirona ORTHOSPHOS XG SL Ready

Phosphor Plate

ScanX (Models)

Classic F3700

Ortho F3750

I/O F3600

Duo D1000F

Swift G8000

Intraoral Cameras

Acteon Sopro (Soprolife, Soprocare, Sopro 717,

Sopro 617)

Air Techniques CamX Polaris and Spectra

intraoral cameras

Kodak/Carestream 1500

intraoral camera

Digital Doc Iris intraoral

USB 2.0 camera

Digital Doc Iris HD USB 3.0

intraoral camera

Schick USB Cam 4 intraoral

camera

Bridges to Other Imaging Software

Apteryx Generic (Other OEM's include; Prof. Suni, Dr. Suni,

Image DDS, Clini View, Clini Touch,

MPDx, VisionDent, Belmont XV,

Belmont XVLite, DentiMax Advance Imaging, DentiMax Imaging,

MaxiViewer, MaxiViewer Lite,

Lightyear, and QSI Image)

CADI (Synca)

CaptureLink

CDR-DICOM (Sirona)

DEXIS

And more

DTX Studio Bridge

Kodak Dental Imaging (Carestream)

Patterson Imaging

Romexis (Planmeca)

RadioVision Bridge

Sidexis (Sirona)

Tigerview (Televere)

Visix (Air Techniques)

Vixwin (Gendex)

XDR (Cyber Medical Imaging)

XrayVision and XVIite (Apteryx Main

imaging software)